

Experiment - 9


Aim :- Spectrophotometric determination of iron with 1, 10, phenanthroline.

Apparatus: Burette, volumetric flasks (50 ml),
cuvettes, funnel, burette, stand
clamp and colorimeter

Chemicals: Mohr's salt salt ($\text{FeSO}_4(\text{NH}_4)_2 \cdot 8\text{O}_4\text{H}_2\text{O}$) 1,10-phenanthroline, acetic acid solution acetate buffer of $\text{PH} 4.5$ and buffer of $\text{PH} 4.5^-$ and H_2SO_4

Indicator : Phenanthroline.

End color: reddish-orange colour



$n = \text{no. of phenanthroline}$

Fe^{2+}

molecules reacting with

Experiment : A spectrophotometric determination of iron
 (II) with 1, 10 - phenanthroline.

Apparatus : Burette, volumetric flasks (50 ml), conical funnel, burette stand, clamp and colorimeter.

Chemicals : Mairis salt solution (ferrous ammonium sulphate $\text{FeSO}_4 (\text{NH}_4)_2 \text{SO}_4 \cdot 6\text{H}_2\text{O}$), 1, 10-phenanthroline, hydroxylamine hydrochloride, acetic acid - sodium acetate buffer of pH 4.5 and sulphuric acid (H_2SO_4)

Theory : Iron (II) reacts with 1, 10 phenanthroline to form an oxime red complex $[\text{C}_{12}\text{H}_{18}\text{N}_2]^+$. The colour intensity is independent of the acidity in pH range 2-9. If iron (III) is present it can be reduced with hydroxylamine hydrochloride. The absorbance is mentioned at a wavelength of 515 nm. The variation of the concentration of a given colored solution, changes the intensity of the transmitted light. The change in light intensity is measured by the instrument called photocolourimeter / colorimeter, where monochromatic light falls on a sample soln, some light is absorbed and the intensity of the transmitted light is decreased. The decrease

Fe concentration in the flask.

K	0.01N
L	$2.0 \times 10^{-5}\text{N}$
M	$4.0 \times 10^{-5}\text{N}$

N	$6.0 \times 10^{-5}\text{N}$
O	$8.0 \times 10^{-5}\text{N}$
P	$10.0 \times 10^{-5}\text{N}$

Observation

- ii) Absorbance of the solution at highest concentration ($10 \times 10^{-5}\text{N}$) at various

S. No	Wavelength (nm)	Absorbance
P ₁	400	
P ₂	420	0.10
P ₃	470	0.19
P ₄	500	0.17
P ₅	530	0.20
P ₆	620	0.14
P ₇	660	0.01
P ₈	706	0.00
		0

$$\lambda_{\text{max}} = 500\text{ nm}$$

The intensity of light is proportional to the thickness of the absorbing medium and the concentration of solution. This may be expressed by Lambert's law.

where c is the concentration of the solution expressed in mol/l and ' E ' is a constant characteristic of the solute and the wave-length of light, E is called the molar extinction coefficient. A is absorbance or optical density (D) of solution, b is path length and is related to the transmittance ($T = I/I_0$)

Procedure -


A - PREPARATION OF SAMPLES -

1. Take six 50 ml volumetric flasks and add 0, 1, 2, 3, 4, 5 ml of FAS solution in each flask. Let's name these volumetric flasks as K, L, M, N, O and P.
2. Then add 2 ml of 1, 10 - phenanthroline solution to each of those volumetric flask.
3. Now dilute each volumetric flask with deionised water to afford a total volume of 50 ml (by filling the these flasks upto the mark) stopper the flasks and mix the contents well by shaking.

(iii) Absorbance of the solution at different concentration at λ_{max}

Sample	Concentration (N)	Absorbance
K	0	0
L	$1 \cdot 7 \times 10^{-5}$	0.04
M	$3 \cdot 7 \times 10^{-5}$	0.09
N	$5 \cdot 7 \times 10^{-5}$	0.13
O	$7 \cdot 7 \times 10^{-5}$	0.17
P	$9 \cdot 7 \times 10^{-5}$	0.20
X	unknown	0.12

A solution with Fe concentration of $2 \times 10^{-5} N$ contains 10 mg of Fe (as used in this exp)



The concentration of unknown sample comes out to be $5 \cdot 2 \times 10^{-5} N$


The Fe content in the unknown given sample is $26 \cdot 1$ mg of Fe ($(e \frac{10}{2 \times 10^{-5}} + 5 \cdot 2 \times 10^{-5})$) mg of Fe

vigorously for few minutes allow the solution to stand for 10 minutes.

4. The first volumetric flask to which ml of FAS is added (i.e. - no Fe^{2+} ions), will serve as a blank (solution K).

B To determine the λ_{\max}

1. Get the two cuvettes issued from the laboratory staff.
2. Fill one of them with the blank solution (K) and another one of the samples containing Fe. Let's say solution P.
3. Light of single wavelength can be produced by selecting the filter on the photochlorimeter. usually the range goes from 410 nm - 700 nm
4. Set the filter to 410nm. Place the cuvette with blank solution, K, in the sample holder.
5. Set the absorbance to 0 %
6. Now place the second cuvette filled with blank solution, K, in the sample holder measure the absorbance of the salt. Now you have the



absorbance of 410 nm for sol. P.

7. By changing the filter to next wavelength each time, repeat steps 4-6. You need to set the absorbance to zero with blank (k) every time you change the + with filter.

8. Now, you have absorbance of solution P, over a range of wavelength from 410 nm - 706 nm you will notice that graph between the absorbance and wavelength takes an inverse parabola shape, with a maximum absorbance around 500 nm or 480 nm. This is λ_{max} .


[Measurement for Absorbance for Solutions L to P at λ_{max} .

9. Set the filter to λ_{max} obtained in part 8.
10. Set the absorbance to 0 using your blank space (k).

Measure the absorbance for solutions L to P now at λ_{max} . Don't disturb the filter in between.

11. Now measure the absorbance for an unknown sample provided to you.

12. Plot absorbance vs concentration for samples L to P.



Connecting maximum points, draw a straight line ideally passing through origin.

6. Using absorbance value for the unknown find out its concentration.

(GRAPH OBSERVATION) : The concentration (x) can be determined from its absorbance value from the graph.

A solution with Fe concentration of 200×10^{-5} M contains 10.06 ug of Fe (as used in this exp) Hence knowing the Fe, concentration of unknown samples, their Fe content in ug can be calculated.

RESULT : The Fe content in the unknown / given sample is 26.1 ug of Fe.